


20.9 Sordaria elliptical spores 25.5 - 28 x 15.5 - 16  
 ascus 22-23 diam. Superba ✓

23.9 Selosa Podospore multispored  
 500 diam. setae lifted, not inflated,  
 but aggregated, < 400 long.  
 transverse body ± immersed  
 sp 14.5 x 9.5-10  
 -16.5  
 reduced st. clavate at tip  
 6 x 2  
 spores in spiral line  
  
 make ascus 190 x 40, with  
 make saccate 260 x 110  
Curvella (H)

27.9 Piptosporium lepidota ✓  
 heterospore  
 dry; glucose hard cells, 2-spored

P. curvella Asubolus primula yellow/green  
 hot seawater on 5.10  
 See 5.11

Thelobolus Stetzneri ✓ (H)

Sordaria fimicola ✓ sp 19.5-20.5 x 10.2-10.5-11  
 ascus 16-18 wide good set

Sord. Superba 25.5-27 x 16

5.10 Podospore asymmetric lifted setae, polysporous. sp 15.5-16 x 9.5  
Curvella ✓

S. limicola ✓

13.10

P. consocella ✓

S. superba ✓ sp

-29  
26 x 16

26.10

5.11

P. consocella ✓

S. superba ✓

Ascobars spores from 200 cultured apothecium



14.5 x 7-8  
-16

crematis ✓

Dried 5.11 with good  
P. consocella + S. superba